



MARINE TESTING SOLUTIONS

Dipslides

Tel: +44 (0)1903 778718
enquiries@marinetestingsolutions.com
www.marinetestingsolutions.com

Marine Testing Solutions Ltd
Unit C1, Rudford Ind Estate
Ford Road
Ford
West Sussex
BN18 0BD



Dipslides

Agar - Dipslides

- Determine aerobic and anaerobic bacteria levels.
- Early indication of bacteria proliferation.
- Results in 24 - 72 hours.
- Inexpensive.
- Easy to use.
- Excellent for trend analysis.

A full range of agar dipslides is available for semi-quantitative determination of aerobic bacteria populations in industrial and recreational waters.

Dipslide accuracy is limited due to the small sample size, but if used correctly and incubated at a constant temperature using the MTS dipslide incubator, they are excellent for the trend analysis and can give an early indication of bacteria proliferation.



- Dipslides are supplied in a cardboard carton containing 10 slides.
- Dipslides have a working shelf life of 6 - 12 months.
- Results are quantified by comparison to a standard density chart supplied.

Description	Product Code
BT2 Aerobic Bacteria	MTS-B100033
SRB Sulphate Reducing Bacteria	MTS-B100034
RBS Aerobic Bacteria	MTS-B100190
SC2 Coliform and E.coli	MTS-B100860
NRB Nitrite Reducing Bacteria	MTS-B100346
TP Aerobic Bacteria / Pseudomonas	MTS-B100032